




FACIAL NERVE DISSECTORS

Stimulate As You Dissect
The perfect warning of nerve proximity

- ✓ Allows surgeon to stimulate while dissecting
- ✓ Suggested for Thyroid and Parotid surgery
- ✓ Removes the need for a separate stimulating probe
- ✓ Optimum tips for precision work
- ✓ Suitable for use with multi-channel nerve monitors
- ✓ Individually packaged and sterile

Neurosign's facial nerve dissectors are designed to stimulate nerves whilst dissecting, removing the need for a separate stimulating probe. Warning of nerve proximity whilst dissecting reduces the risk of facial nerve damage and improves efficiency by reducing surgical time spent changing instruments.




Contact your Neurosign® representative for further information.

www.neurosignsurgical.com
info@neurosignsurgical.com



REDUCING THE RISK OF FACIAL NERVE DAMAGE WHILST IMPROVING SURGICAL EFFICIENCY

- ✓ **Reduces Risk** - stimulating whilst dissecting reduces the risk of facial nerve damage
- ✓ **Improves Efficiency** - removes the need for a separate stimulating probe reducing time spent changing implements.
- ✓ **Cost effective** - single use instrument provided individually packaged and sterile
- ✓ **Compatible** - suitable for use with multi-channel nerve monitors



FACIAL NERVE DISSECTOR ORDERING INFORMATION:

4013-00

Facial Nerve Dissectors, box of 10,
Sterile, single use, individually packaged.

NMK1709-01

© 2017, The Magstim Company Limited. Magstim® and Neurosign® are registered trademarks of The Magstim Company Limited. All specifications are subject to change. The information contained within this document is accurate at the time of publication.

 Technomed
Amerikalaan 71, 6199 AE Maastricht-Airport, The Netherlands

 **NEUROSIGN**®
www.neurosignsurgical.com